

Mist Shower Unit

The Mist shower is a self-contained chamber that sprays a fine mist inside the chamber. ALLEGRO manufacture Mist Showers to ensure that gowns are wet from all the sides. We incorporate intelligent controls to ensure the exact / correct time is allocated to wet the dress.



Technical Specification:

- Air Cleanliness - ISO 5 / Class 100
- Solid Door with Electromagnet Interlocking
- SS 304, Mist, Water and Air nozzles to ensure complete coverage of inside
- Electrical solenoid valves for water and air lines used for mist generation
- Door panels interlocked with emergency push button for exit
- HEPA / ULPA filtration
- PAO test port at the upstream side of the HEPA filter
- Union Connector for Oxygen Input
- Perforated floor panel to drain