

Mobile LAF (class - 100 trolley)

The self-powered unit by means of dry UPS battery provides a laminar air flow of class ISO 5. ALLEGRO Integrate the Mobile LAF with the upstream and downstream equipment for easy of transferring the material.



Technical Specification:

- Cabinet fabricated in SS 304/316
- Air flow direction: Vertical / Horizontal
- Design according to GMP-guideline easy to clean and to disinfect.
- PU wheels with brakes for transportation
- Specially design centrifugal or back-ward curved type air blower
- Magnehelic differential pressure gauge
- Soft touch LED controller for blower/light/pressure display
- In build hour meter
- Mini-plate HEPA filters H-14, 99.997% down to 0.3μ.
- Pre-filter non-woven, HDPE mesh EU-04 95% efficiency 5 μ
- UV light.
- Florescent light
- PAO port/atmospheric nozzles
- Sound level minimum 65db on scale
- Pre-Installed battery.
- Power supply single phase 220v 50 HZ

Model no.	HEPA Size L x W	Over all Dimension W x D x H
AGPL/MLFR01	610 x 610	760 x1170 x2240
AGPL/MLFR02	610 x 915	1070 x1170 x 2240
AGPL/MLFR03	610 x 1220	750 x 1370 x 1950