

Static Pass Box

Especially design Static Pass Box for transfer of material and to maintain class between to areas. These static pass boxes are installed in manufacturing areas and filter the incoming air in order to keep the air free from dust and other contaminants with UV light along with hour metre.



Technical Specification:

- Cabinet fabricated in SS 304/316
- Design according to GMP-guideline easy to clean and to disinfect.
- Interlocking door system
- Soft touch LED controller for display
- In build hour meter
- UV light.
- Florescent light
- Power supply single phase 220v 50 HZ

Model no.	Working Area W x D x H	Over all Dimension W x D x H
AGPL/SPB01	200x250x200	305x305x305
AGPL/SPB02	350x400x350	460x460x460
AGPL/SPB03	500x550x500	610x610x610
AGPL/SPB04	650x710x650	760x760x760