

Sterile Garment Cubicle

ALLEGRO is a superior manufacturer of Sterile Garments Cubicle used as storing cubicles or cabinets which are used in pharmaceutical cleanrooms for operator gowning and sterile cloths switch storage material and aids in preserving a sense of class between the two regions. It is used between the two different classified areas.



Technical Specification:

- Cabinet fabricated in SS 304/316
- Air flow direction: Vertical
- Design according to GMP-guideline easy to clean and to disinfect.
- Interlocking door system
- Specially design centrifugal or back-ward curved type air blower
- Magnehelic differential pressure gauge
- Soft touch LED controller for blower/light/pressure display
- In build hour meter
- Mini-plate HEPA filters H-14, 99.997% down to 0.3μ.
- Pre-filter non-woven, HDPE mesh EU-04 95% efficiency 5 μ
- UV light.
- Florescent light
- PAO port/atmospheric nozzles
- Sound level minimum 65db on scale
- Power supply single phase 220v 50 HZ

Model no.	Working Area W x D x H	Over all Dimension W x D x H
AGPL/SGC01	900 x 450 x 1220	915 x 460 x 1680 (wall mount)
AGPL/SGC02	900 x 450 x 1830	915 x 460 x 2400 (floor mount)